



BROWN

DECONTAMINATING EQUIPMENT

To be completed by lab personnel:

Check One: Disposal On Campus Relocation Repair Storage

Equipment:

Location (Building and Room Number):

Laboratory Supervisor and telephone number:

Hazard Reduction Checklist:

Biological material disinfected/decontaminated with effective agent listed in the Exposure Control Plan (ECP): 10% Bleach Soap and Water Other: _____

Biosafety Cabinets must be disinfected by a contracted company like B & V Testing 1-800-851-9081.

Equipment was never used for radioactive materials.

If used with radioactive materials, the equipment has been fully monitored for radioactive contamination and has been decontaminated. Forward a copy of the wipe test results to Box 1914. Attn. RSO.

Cleared of all other recognizable hazards (broken glass, razors, other sharps or medical waste).

No knowledge of PCB oils, mercury switches, thermometers or other hazardous materials within.

All hazard labels removed or properly defaced upon completion of decontamination of equipment.

Signature of personnel that performed above procedures

Date

Once all of the above procedures have been checked off, post a copy directly on the piece of equipment to communicate to workers or movers the equipment is safe to handle. When equipment is being disposed also include another copy of this form with the Brown University Surplus Disposition Form. For questions regarding this form please contact Environmental Health & Safety at 863-3353.

**THIS IS NOT A SURPLUS RELEASE FORM. ALL SURPLUS EQUIPMENT MUST BE REPORTED TO:
Surplus@Brown.edu**